Electra S4

ANTI-SPASM

I-TECH
MEDICAL DIVISION

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Product description

1. Remote control
2. ON/OFF Switch
3. Pedals (2 pcs - disassembled at the time of purchase)
4. Handles (2 pcs - disassembled at the time of purchase)
5. Bases (2 pcs – disassembled at the time of purchase)
6. Component mounting key

Operational characteristics

- Anti-spasm
- Remote control with on-board computer
- Adjustable speed: 15-30-45-60 revs/minute
- Setting time: 3-6-9-12-15 min
- Forward thrust on pedals
- Reverse thrust on pedals
- Fixing of pedals and handles made easy
- Adjustable thrust on pedal radius: 13-11-9 cm
- Overload motor protection with automatic stop
- Tachometer function
- Calories counter function
- Countdown function
- Pedals or handles operation
- Preset training function (auto)
- Data and functions readable on the electronic display
- Pierced base supports for vertical position of appliance
- Top handle bar

Anti-spasm and user safety

Electra S4 is equipped with the new muscular spasm survey system, design by I.A.C.E.R. Medical Division R&D Department to guarantee user safety and an efficient and safe training.

During the movement a quick and unexpected muscular spasm can occur, and the muscles will rapidly tense.

Electra S4 recognizes the spasm and stops pedals/handles rotation for 3 seconds, and then inverts the rotation to untie muscular spasm.

The spasm is untied following the antagonist inhibition principle: a bending spasm is untied through muscular extension, while an extension spasm is untied through bending.

If a second spasm will occur, Electra S4 stops and warns the user with an acoustic signal: in this way the user safety is guaranteed, in particular for tetraplegic patients who need to be followed by a familiar or a doctor. To restart Electra S4 simply press the start/pause button.

Technical characteristics

Power supply. 230 V - 50 Hz. **Electrical motor.** 230 V - 50 Hz.

Motor revolutions. 48 min. **Motor power.** 50 W.

Timer. 3 - 15 min. **Pedal semi diameter.** 13 - 11 - 9 cm.

Weight. 5.7 kg. **Dimensions.** 440 x 425 x 285 mm.

Adjustable thrust on pedal radius (fig. 3). Thanks to the new support, pierced in three different positions, pedals or handles can be mounted by choosing the rotation radius more suitable with the height of the user and therefore adapt the movement extension to the mobility of the joints. The radii are as follows: external hole R=13 cm, intermediate hole R=11 cm, internal hole R=9 cm. The radius variations apply to the subjects, respectively: Tall (R=13), Medium (R=11) and Short (R=9).

Thrust on pedal speed (fig. 5-6-7). Four rotation groups suitable for any rehabilitative requirement: Low = 15 RPM, Medium-low = 30 RPM, Medium-high = 45 RPM, High = 60 RPM.

Fixing of pedals and handles made easy (fig. 1). Insert the pivot pin of the pedal or of the handle in the most suitable hole for the user's height and lock it by screwing the supplied bolt.

Instructions for using the appliance

- A. Choose if you wish to use the pedals (fig. 2)
- B. Choose if you wish to use the handles (fig. 4)
- C. Insert the appliance plug in a 230/50 socket
- D. Start the appliance by using the ON/OFF switch
- E. When you start the appliance, four lines of symbols and numbers will appear on the remote control display (fig. 5-6-7) which allow to display and select the training programmes
- F. Press the START/PAUSE button to start, interrupt and restart the training
- G. At the end of the training, switch the appliance off by using the ON/OFF switch

Reading the display

FIGURE 5

AUTO. Preset training function activated and deactivated through AUTO button. A line below the symbol shows that the function is on.

15 30 45 60 RPM. Rotation revs per minute. A line below the number shows the selected speed.

3 6 9 12 15 MIN. Training time which can be set in minutes. A line below the number shows the selected time.

0 COUNT. Shows the number of revs performed by the machine. To select the following function (fig. 6) press the MODE button.

NOTE: it is possible to modify all the parameters while the machine is running.

FIGURE 6

0.0 CAL. Shows the CALORIES burnt by training. To select the following function (fig. 7) press the MODE button.

FIGURE 7

TIME. Shows the active exercise time in the machine. To modify it, press the TIMER button scrolling the line up to below the number.

WARNINGS

Persons not in full possession of their faculties **are not allowed** from using the machine. Even temporarily disabled people **are not allowed** from using the machine unless they are assisted by qualified personnel (ex. Doctor or Therapist). It is **not allowed** to use the appliance if it shows wear and tear signs such as damaging of the plug or power cable after violent falls or if it becomes damp by contact with any liquid. It is **not allowed** to use the appliance near inflammable substances or in environments with high oxygen concentrations. It is **not allowed** to let the appliance connected to the power source when it is not in use or during the replacement of Pedals and/or Handles. It is **not allowed** to use the appliance below sheets or pillows; overheating could provoke fire and/or short-circuits causing injuries to persons. It is **not allowed** to use the appliance if ventilation slots are obstructed. It is **not allowed** to use the appliance outdoor.

Safety note


The appliance is equipped with a protection against incorrect use: in case of overloading or incorrect use, the motor automatically stops. An acoustic signal warns the user. To restart the appliance, press the ON/OFF switch. If this operation is not carried out, after one minute the acoustic signal automatically stops.

Maintenance

Cleaning. When required, clean the appliance with a dry cloth. Never dip the appliance into the water or use solvents or chemical products.

Maintenance. Never roll the power cable around the appliance and keep it in a safe and dry place.

Remote control

1. Start/Pause. Start, stop and pause button
2. Display
3. Timer. Presetting training time button
4. Speed. Rotation speed adjustment button
5. Auto. Preset training selection button
6. Forward/Reverse. Button for direction of rotation selection
7. Mode. Data display button

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IACER Srl